

Urology

Interventional

Biopsy Solutions

Breast care

Radiology


Tissue Marker

R&D

vigeo

Export Catalogue

Medical Device Regulation  
EU 2017/745



vigeo



**SOFT TISSUE BIOPSY**

COMPACT	<i>pag.4</i>
COMPACT T	<i>pag.5</i>
COMPACT URO	<i>pag.6</i>
COREX	<i>pag.7</i>
VTEK	<i>pag.8</i>
VVTEK	<i>pag.9</i>
AUTOCUT	<i>pag.10</i>
CT-CORE	<i>pag.11</i>
INTRO/INTRO PLUS	<i>pag.12</i>
HEPASHOT	<i>pag.13</i>

---

**SOFT TISSUE MARKING**

BIOMARC	<i>pag.14</i>
BIOMARC IGRT	<i>pag.15</i>

---

**WOMEN'S HEALTH**

RELOCK PREMIUM SOFT	<i>pag.16</i>
RELOCK	<i>pag.17</i>
RELOCK DUO	<i>pag.18</i>
RELOCK PREMIUM	<i>pag.19</i>
RELOCK PREMIUM DUO	<i>pag.20</i>
SURGIMARC	<i>pag.21</i>
SURGIMARC SLIM	<i>pag.22</i>
SCREEN VISION	<i>pag.23</i>
ISTERIX	<i>pag.24</i>

---

**DRAINAGE**

V-DRAIN	<i>pag.25</i>
---------	---------------

**FINE NEEDLE ASPIRATION**

CYTOCAN	<i>pag.26</i>
PRENATAL	<i>pag.27</i>

---

**CLINICAL ONCOLOGY**

BLOCK	<i>pag.28</i>
BCUT	<i>pag.29</i>
MYELO SLIM	<i>pag.30</i>
MYELO	<i>pag.31</i>
ERGO	<i>pag.32</i>

---

**SKIN CARE**

D PUNCH/CURETTE	<i>pag.33</i>
-----------------	---------------

---

**INFECTION CONTROL**

DECOMAT	<i>pag.34</i>
---------	---------------

# COMPACT SOFT TISSUE BIOPSY


Automatic soft tissue biopsy device

CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )



- One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer
- Programmable: releases the stylet allowing time for tissue to settle in notch (ideal for denser tissue or when there is a desire to verify notch position) and activates the cutting cannula
- Thin-wall cannula for larger core
- Colour coded for quick gauge size identification
- Depth markings etched on cannula for precise positioning
- The Coaxial kit permits multiple biopsies (US and CT guided)
- Loaded status indicator
- Plastic safety protection to avoid accidental firing



Colour Code

One-handed design with loaded status indicator

## COMPACT Qty/Box 10pcs

CODE	SIZE	COAXIAL KIT
● VCT1210	12Gx10cm	VCTC1210
● VCT1410	14Gx10cm	VCTC1410
● VCT1411	14Gx11cm	VCTC1411
● VCT1413	14Gx13cm	VCTC1413
● VCT1415	14Gx15cm	VCTC1415
● VCT1610	16Gx10cm	VCTC1610
● VCT1615	16Gx15cm	VCTC1615
● VCT1620	16Gx20cm	VCTC1620
● VCT1810	18Gx10cm	VCTC1810
● VCT1815	18Gx15cm	VCTC1815
● VCT1820	18Gx20cm	VCTC1820
● VCT1825	18Gx25cm	VCTC1825
● VCT2010	20Gx10cm	VCTC2010
● VCT2015	20Gx15cm	VCTC2015

On request Vigeo is available to provide needles with different diameter and length

# COMPACT T SOFT TISSUE BIOPSY


Short notch automatic soft tissue biopsy device

CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )



- One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer
- Programmable: releases the stylet allowing time for tissue to settle in notch (ideal for denser tissue or when there is a desire to verify notch position) and activates the cutting cannula
- Sample notch 10mm for critical procedures
- Thin-wall cannula for larger core
- Colour coded for quick gauge size identification
- Depth markings etched on cannula for precise positioning
- Loaded status indicator
- Plastic safety protection to avoid accidental firing



Colour code and plastic safety protection

## COMPACT T Qty/Box 10pcs

CODE	SIZE
● VCT11410	14Gx10cm
● VCT11413	14Gx13cm
● VCT11415	14Gx15cm
● VCT11810	18Gx10cm
● VCT12010	20Gx10cm


On request Vigeo is available to provide needles with different diameter and length

# COMPACT URO SOFT TISSUE BIOPSY

Automatic soft tissue biopsy device

CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )



- One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer
- Thin-wall cannula for larger core
- Colour coded for quick gauge size identification
- Depth markings etched on cannula for precise positioning
- Loaded status indicator
- Plastic safety protection to avoid accidental firing



Flat design for transrectal ultrasound (TRUS)-guided needle biopsy



Colour code

## COMPACT URO Qty/Box 10pcs

CODE	SIZE
● VCTA1620	16Gx20cm
● VCTA1625	16Gx25cm
● VCTA1820	18Gx20cm
● VCTA1825	18Gx25cm
● VCTA1830	18Gx30cm

On request Vigeo is available to provide needles with different diameter and length





- One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer
- Equipped with a visual and mechanical safety-system
- Sample notch of 19 mm
- Thin-wall cannula for larger core
- Two triggers permit ergonomic handling for right and left-handed application, regardless of the position of the extraction location
- Colour coded for quick gauge size identification
- Depth markings etched on cannula for precise positioning
- The Coaxial kit permits multiple biopsies (US and CT guided)



Mechanical safety-system



Longer handle for ergonomicity improvement

### COREX Qty/Box 10pcs

CODE	SIZE	COAXIAL KIT
● VBS1410	14Gx10cm	VBSC1410
● VBS1415	14Gx15cm	VBSC1415
● VBS1610	16Gx10cm	VBSC1610
● VBS1615	16Gx15cm	VBSC1615
● VBS1620	16Gx20cm	VBSC1620
● VBS1625	16Gx25cm	VBSC1625
● VBS1810	18Gx10cm	VBSC1810
● VBS1815	18Gx15cm	VBSC1815
● VBS1820	18Gx20cm	VBSC1820
● VBS1825	18Gx25cm	VBSC1825

On request Vigeo is available to provide needles with different diameter and length



- Safety button to prevent unwanted shots
- Two triggers
- Two possible puncture depths of 16 or 22 mm
- Compatible exclusively with Vigeo Autocut biopsy needle
- Loaded status indicator



Loaded status indicator



Coaxial kit ( VGC )

VTEK		
CODE	THROW	SAFETY
VTEK	16-22 mm	YES

AUTOCUT Qty/Box 10pcs		
CODE	SIZE	COAXIAL KIT
● VVG1210	12Gx10cm	VGVC1210
● VVG1410	14Gx10cm	VGVC1410
● VVG1413	14Gx13cm	VGVC1413
● VVG1415	14Gx15cm	VGVC1415
● VVG1610	16Gx10cm	VGVC1610
● VVG1615	16Gx15cm	VGVC1615
● VVG1620	16Gx20cm	VGVC1620
● VVG1810	18Gx10cm	VGVC1810
● VVG1815	18Gx15cm	VGVC1815
● VVG1820	18Gx20cm	VGVC1820
● VVG1825	18Gx25cm	VGVC1825
● VVG2010	20Gx10cm	VGVC2010
● VVG2015	20Gx15cm	VGVC2015

On request Vigeo is available to provide needles with different diameter and length





- One trigger
- Puncture depth of 22 mm
- Compatible exclusively with Vigeo Autocut biopsy needle
- Loaded status indicator



Loaded status indicator



Coaxial kit ( VGCC )

VVTEK		
CODE	THROW	SAFETY
VVTEK	22 mm	NO

AUTOCUT <small>Qty/Box 10pcs</small>		
CODE	SIZE	COAXIAL KIT
● VVG1210	12Gx10cm	VGCC1210
● VVG1410	14Gx10cm	VGCC1410
● VVG1413	14Gx13cm	VGCC1413
● VVG1415	14Gx15cm	VGCC1415
● VVG1610	16Gx10cm	VGCC1610
● VVG1615	16Gx15cm	VGCC1615
● VVG1620	16Gx20cm	VGCC1620
● VVG1810	18Gx10cm	VGCC1810
● VVG1815	18Gx15cm	VGCC1815
● VVG1820	18Gx20cm	VGCC1820
● VVG1825	18Gx25cm	VGCC1825
● VVG2010	20Gx10cm	VGCC2010
● VVG2015	20Gx15cm	VGCC2015

On request Vigeo is available to provide needles with different diameter and length

Refill needle for reusable biopsy system



**AUTOCUT V (VVG)** Specially designed to use with VTEK - VVTEK Vigeo system - MAGNUM® Bard system



**AUTOCUT P (VVG P)** Specially designed to use with VVTEK Vigeo system - PRO-MAG® ULTRA Argon system



**AUTOCUT B (VVE)** Specially designed to use with PRO-MAG® 2.2 Argon system - BIOPTY® Bard system

- Single patient core biopsy needle designed for use with reusable biopsy system
- Available in several needle gauge sizes and lengths
- Centimeter depth markings assist in needle placement
- Attached removable spacer (VVG / VVG P) maintains the distance between the two needle hubs promoting ease of handling and simplifying targeting
- The coaxial outer needle permits multiple biopsies (US and CT guided)
- Colour coded for quick gauge size identification

### AUTOCUT Qty/Box 10 pcs

Add **C** for Coaxial kit ( VVGC / VVGPC / VVEC )

Add **1** for Trocar tip ( VVG1 / VVG P1 )

CODE	SIZE
● VVG1210	12Gx10cm
● VVG1410 VVG P1410 VVE1410	14Gx10cm
● VVG1413 VVG P1413 VVE1413	14Gx13cm
● VVG1415 VVG P1415 VVE1415	14Gx15cm
● VVG1610 VVG P1610 VVE1610	16Gx10cm
● VVG1615 VVG P1615 VVE1615	16Gx15cm
● VVG1620 VVG P1620 VVE1620	16Gx20cm
● VVG1810 VVG P1810 VVE1810	18Gx10cm
● VVG1815 VVG P1815 VVE1815	18Gx15cm
● VVG1820 VVG P1820 VVE1820	18Gx20cm
● VVG1825 VVG P1825 VVE1825	18Gx25cm
● VVG2010 VVG P2010 VVE2010	20Gx10cm
● VVG2015 VVG P2015 VVE2015	20Gx15cm

On request Vigeo is available to provide needles with different diameter and length

Semiautomatic guillotine biopsy device



- Lightweight and compact to facilitate one-handed operation during procedures
- Precise, rapid-firing mechanism for an intact core tissue sample
- Specimen notch fully exposed when stylet is advanced for ease of sample removal
- The puncture depth can be adjusted between 15 to 22 mm, the notch length between 10 to 20 mm
- Beveled stylet for easy penetration into the lesion with less trauma to surrounding tissues and organs
- Colour coded for quick gauge size identification
- Echogenic tip for enhanced visualization under ultrasound guidance
- Depth markings etched on cannula for precise positioning
- The Coaxial outer needle permits multiple biopsies (US and CT guided)



One-handed design



Coaxial kit ( VVD0C )

### CT-CORE Qty/Box 10pcs

CODE	SIZE	COAXIAL KIT
● VVD01410	14Gx10cm	VVD0C1410
● VVD01415	14Gx15cm	VVD0C1415
● VVD01610	16Gx10cm	VVD0C1610
● VVD01615	16Gx15cm	VVD0C1615
● VVD01620	16Gx20cm	VVD0C1620
● VVD01810	18Gx10cm	VVD0C1810
● VVD01815	18Gx15cm	VVD0C1815
● VVD01820	18Gx20cm	VVD0C1820
● VVD01825	18Gx25cm	VVD0C1825
● VVD02010	20Gx10cm	VVD0C2010
● VVD02015	20Gx15cm	VVD0C2015

On request Vigeo is available to provide needles with different diameter and length

# INTRO/INTRO PLUS SOFT TISSUE BIOPSY

CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )

Standard coaxial introducer needle used for single access in multiple biopsies



- Echogenic treatment for ultrasound guided procedures
- Centimeter marking and optional sliding stopper
- Sterilized by ETO, shelf life 5 years
- INTRO (code VVT) with blunt cannula and trocar tip
- INTRO PLUS (code VVTC) with sharp cutting cannula and trocar tip



Depth marker with code colour

## INTRO/INTRO PLUS Qty/Box 20pcs

CODE	SIZE
● VVT1105 VVTC1105	11Gx5cm
● VVT1305 VVTC1305	13Gx5cm
● VVT1310 VVTC1310	13Gx10cm
● VVT1315 VVTC1315	13Gx15cm
● VVT1505 VVTC1505	15Gx5cm
● VVT1510 VVTC1510	15Gx10cm
● VVT1515 VVTC1515	15Gx15cm
● VVT1705 VVTC1705	17Gx5cm
● VVT1710 VVTC1710	17Gx10cm
● VVT1715 VVTC1715	17Gx15cm

On request Vigeo is available to provide needles with different diameter and length

One-handed Menghini aspiration device



- One-handed safe and fast operation
- Pyramidal tip allows for easy insertion
- Menghini design cannula for cylindrical and large sampling
- Intended for use in liver, lymph nodes and thyroid
- Colour coded for quick gauge size identification
- Centimeter depth markings assist in needle placement



One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer

## HEPASHOT Qty/Box 10pcs

CODE	SIZE
● VVC01610	16Gx10cm
● VVC01615	16Gx15cm
● VVC01620	16Gx20cm
● VVC01805	18Gx5cm
● VVC01810	18Gx10cm
● VVC01815	18Gx15cm
● VVC01820	18Gx20cm
● VVC02005	20Gx5cm
● VVC02010	20Gx10cm
● VVC02015	20Gx15cm
● VVC02020	20Gx20cm
● VVC02105	21Gx5cm
● VVC02110	21Gx10cm
● VVC02115	21Gx15cm

On request Vigeo is available to provide needles with different diameter and length


## Breast tissue marker device



- Barbell shaped pyrolytic carbon coated ceramic with exceptionally biocompatibility and clear imaging characteristics
- Unique shapes designed to prevent migration
- No allergy response or hypersensitivity reported
- Designed for safe, reliable and precise placement during core biopsies
- Permanent and clear visibility with mammography, ultrasound and MRI
- No interference with MRI visualization, lower level of artifact than metal markers
- No interference with tomosynthesis visualization
- Available in three distinct size configurations
- Available for vacuum assisted and core biopsy needle devices

### WHEN TO USE

- Marking after removal of tumor
- Marking microcalcification area after breast vacuum assisted biopsy (VAB)
- Marking of suspicious breast lesions after critical core biopsy
- Orientation for planning radiation therapy



Ultrasound



Mammography

BIOMARC Qty/Box 10pcs			
CODE	SIZE	MARKER	DESCRIPTION
040220	8Gx165,3mm	2x4mm	Flex. Applicator, 8G Mammotome®
040203	9Gx165,3mm	2x4mm	Flex. Applicator, 11G Mammotome® / 9G Eviva®
040204	9 G x 1 2 0 m m	2x4mm	Flex. Applicator, 9G Eviva® / 10G Spirotome®
040239	9Gx363,7mm	2x4mm	Flex. Applicator, 10G SenoRX® / 10G Bexcore® / 9GAtec®
040306	13Gx117,7mm	2x4mm	Rigid Applicator, 13G
040320	14Gx117,7mm	1x5mm	Rigid Applicator, 14G
040300	14Gx117,7mm	1x3mm	Rigid Applicator, 14G
040310	17Gx117,7mm	1x5mm	Rigid Applicator, 17G

On request Vigeo is available to provide needles with different diameter and length




# BIOMARC IGRT SOFT TISSUE MARKING

Exclusive Importer **CE**  
2797




Fiducial marker for radiation therapy




- Used to target localization for lesions or organs which move in respect to external anatomy
- Pure gold cylinders version specially knurled to inhibit migration
- Pyrolytic carbon coated zirconium oxide or tantalum barbell version to reduce the artifacts
- Provide visualization under fluoroscopy, x-ray, CT, MRI, ultrasound and other imaging modalities
- Preloaded in a pre-waxed tipped needle or in a wax-free special version
- Reduces the artifact in the MV setting (BiomarC EX)
- Ideal for proton or photon therapy (BiomarC Pro)

BiomarC GOLD




BiomarC KV ENH



BiomarC KV




BiomarC EX



BiomarC PRO



BiomarC SECURE



## BIOMARC IGRT

CODE	SIZE	DESCRIPTION	Qty/Box
040116	17Gx20cm	1,2x3mm Gold Marker	3
040118	17Gx20cm	1,2x4,2mm Gold Marker-Ergonomic handle	10
040141	18Gx20cm	0,9x3mm Gold Marker	3
040125	17Gx30cm	1,2x3mm Gold Marker	3
040117	18Gx20cm	1x5mm Carbon KV ENH Marker	3
040113	18Gx20cm	1x3mm Carbon KV Marker	3
040146	18Gx20cm	1x5mm Carbon EX Marker	3
040148	18Gx20cm	1x10mm Carbon PRO Marker	1
040150	Sutured	1,7x4,6mm Carbon Secure Marker	10

On request Vigeo is available to provide needles with different diameter and length

# RELOCK PREMIUM SOFT WOMEN'S HEALTH

Preoperative repositionable breast localization device




In compliance with MDR  
( Regulation EU 2017/745 )

PATENT PENDING




- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Retractable hook with unique flexibility of its polyamide guiding thread
- Indicator on wing shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- 9mm loop hook
- Single packaged stiffener cannula used in the operating room




1




1a




2



2a



3



Mammography

## RELOCK PREMIUM SOFT Qty/Box 10pcs


CODE	SIZE
● VVRPS1905	19Gx5cm
● VVRPS1907	19Gx7cm
● VVRPS1910	19Gx10cm
● VVRPS2005	20Gx5cm
● VVRPS2007	20Gx7cm
● VVRPS2010	20Gx10cm

On request Vigeo is available to provide needles with different diameter and length


Preoperative repositionable breast localization device



- Preoperative marking of local lesions
- Retractable hook with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Wire fixation clamp to avoid migration
- 9mm loop hook



Mammography



Wire fixation clamp

### RELOCK Qty/Box 10pcs

CODE	SIZE
● VVR1905	19Gx5cm
● VVR1907	19Gx7cm
● VVR1910	19Gx10cm
● VVR2007	20Gx7cm
● VVR2010	20Gx10cm
● VVR2015	20Gx15cm

On request Vigeo is available to provide needles with different diameter and length

# RELOCK DUO WOMEN'S HEALTH


Preoperative repositionable breast localization device




In compliance with MDR  
( Regulation EU 2017/745 )




- Preoperative marking of local lesions
- Retractable hooks with unique flexibility
- Echogenic needle tip for ultrasound use
- Hooks release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Fixation piece for preoperative fixing of the wire
- Twisted wire for enhanced visibility and lower risk of accidental cutting
- Double 5mm loop hook




Half exposure



Full exposure



Wire fixation clamp



Mammography

## RELOCK DUO Qty/Box 10pcs

CODE	SIZE
VVNSP1910	19Gx10cm
VVNSP2005	20Gx5cm
VVNSP2007	20Gx7cm
VVNSP2010	20Gx10cm
VVNSP2012	20Gx12cm
VVNSP2015	20Gx15cm

On request Vigeo is available to provide needles with different diameter and length

# RELOCK PREMIUM WOMEN'S HEALTH

Preoperative repositionable breast localization device




CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )


PATENT PENDING




- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Removable handle for mammography procedure
- Retractable hook with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- 9mm loop hook



Mammography



Wire fixation clamp  
with protection tube



Removable handle

## RELOCK PREMIUM Qty/Box 10pcs

CODE	SIZE
● VVRP1905	19Gx5cm
● VVRP1907	19Gx7cm
● VVRP1909	19Gx9cm
● VVRP1912	19Gx12cm

On request Vigeo is available to provide needles  
with different diameter and length



# RELOCK PREMIUM DUO WOMEN'S HEALTH

Preoperative repositionable breast localization device


CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )

PATENT PENDING




- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Removable handle for mammography procedure
- Retractable hooks with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hooks release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- Double 5mm loop hook



Half exposure




Full exposure



Mammography



Removable handle



Wire fixation clamp  
with protection tube


## RELOCK PREMIUM DUO Qty/Box 10pcs

CODE	SIZE
VVNSPX1807	18Gx7cm
VVNSPX1810	18Gx10cm
VVNSPX2005	19,5Gx5cm
VVNSPX2009	19,5Gx9cm
VVNSPX2012	19,5Gx12cm

On request Vigeo is available to provide needles with different diameter and length




Preoperative breast localization device



- Preoperative marking of local lesions
- Puncture cannula with ultrasonic tip with guidance wire in place
- Transparent cannula connector
- The Z-shaped hook provides safer delivery in critical areas
- The X-shaped hook provides excellent fixation and is intended to limit migration



Single Hook



Z-shaped hook



X-shaped hook

SURGIMARC Qty/Box 20 pcs			
CODE			SIZE
● VVN1807	VVN1807	VVNZ1807	18Gx7cm
● VVN1810	VVN1810	VVNZ1810	18Gx10cm
● VVN1815	VVN1815	VVNZ1815	18Gx15cm
● VVN2005	VVN2005	VVNZ2005	20Gx5cm
● VVN2007	VVN2007	VVNZ2007	20Gx7cm
● VVN2009	VVN2009	VVNZ2009	20Gx9cm
● VVN2010	VVN2010	VVNZ2010	20Gx10cm
● VVN2015	VVN2015	VVNZ2015	20Gx15cm
● VVN2107	VVN2107	VVNZ2107	21Gx07cm
● VVN2110	VVN2110	VVNZ2110	21Gx10cm

On request Vigeo is available to provide needles with different diameter and length


One-handed preoperative breast localization device



- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Puncture cannula with ultrasonic tip with guidance wire in place
- The Z-shaped hook provides safer delivery in critical areas (available in 18G version)
- The X-shaped hook provides excellent fixation and is intended to limit migration (available in 20 and 21G version)



Single Hook



Z-shaped hook



X-shaped hook




One-handed design

### SURGIMARC SLIM Qty/Box 10pcs

CODE	SIZE
● VVNS1807 VVNSX1807 (Z-shaped hook)	18Gx7cm
● VVNS1810 VVNSX1810 (Z-shaped hook)	18Gx10cm
● VVNS2007 VVNSX2007 (X-shaped hook)	20Gx7cm
● VVNS2010 VVNSX2010 (X-shaped hook)	20Gx10cm
● VVNS2015 VVNSX2015 (X-shaped hook)	20Gx15cm
● VVNS2107 VVNSX2107 (X-shaped hook)	21Gx7cm
● VVNS2110 VVNSX2110 (X-shaped hook)	21Gx10cm


On request Vigeo is available to provide needles with different diameter and length

Tissue transfer and x-ray system



- Safely image and transport of surgical breast specimens without losing orientation
- Extremely radiolucent and radiopaque grid with pinpointing side indicators
- Precise anatomical assignment of sample margins
- Special stands allowing horizontal and vertical examination
- Once pinned the tumour tissue remains stable and can be easily transported
- On request SCREEN VISION SUPPORT for vertical radiographs

AREOLAR   ESTERNAL   DEEP   MEDIAL   CAUDAL   CRANIAL



VVX3P orientation pins version

SCREEN VISION Qty/Box 5pcs	
CODE	THICKNESS
VVX2	20mm
VVX3	30mm
VVX3P	30mm + Pins

On request Vigeo is available to provide needles with different diameter and length

## Hysterosalpingography and sonohysterography catheter



- Flexible, radiopaque polyurethane catheter shaft
- Soft, latex-free concentric balloon
- Atraumatic catheter tip
- Front delivery catheter tip version allows guide wire introduction
- Shapeable version with an embedded stylet in cases when stenosis is present or when the vaginal tract and cervical canal are not aligned
- Stable placement sheath for easy introduction of the catheter into the cervical canal
- Syringe included, preset and mounted on the stopcock for balloon inflation port



Shapeable



Front delivery

### ISTERIX Qty/Box 10 pcs

CODE	SIZE	DESCRIPTION
VMI1154	5F	Side delivery
VMI1155	7F	Side delivery
VMI1185	5F	Shapeable
VMI1188	5F	Front delivery

On request Vigeo is available to provide needles with different diameter and length


Single step drainage catheter




- Large oval drainage side holes
- Large lumen and tapered tip
- Locking system available
- Internal diameter of hub and catheter match
- Centimeter depth markings assist in catheter placement
- Echogenic needle tip with trocar stylet
- Hydrophilic coating for easy penetration
- Easy-to-remove pigtail



Large oval holes



Trocar stylet



Locking system available

### V-DRAIN with locking pigtail Qty/Box 5pcs

CODE	SIZE
VVVB00620	6Fx19cm
VVVB00630	6Fx29cm
VVVB00720	7Fx19cm
VVVB00730	7Fx29cm
VVVB00820	8Fx19cm
VVVB00830	8Fx29cm
VVVB01020	10Fx19cm
VVVB01030	10Fx29cm
VVVB01430	14Fx29cm
VVVB01630	16Fx29cm

### V-DRAIN without locking pigtail Qty/Box 5pcs

CODE	SIZE
VVVB0620	6Fx19cm
VVVB0630	6Fx29cm
VVVB0720	7Fx19cm
VVVB0730	7Fx29cm
VVVB0820	8Fx19cm
VVVB0830	8Fx29cm
VVVB1020	10Fx19cm
VVVB1030	10Fx29cm
VVVB1430	14Fx29cm
VVVB1630	16Fx29cm

On request Vigeo is available to provide needles with different diameter and length

Chiba needle for aspiration biopsy



- Puncture cannula developed especially for ultrasonic diagnostics
- Echogenic needle tip assists precise needle tip positioning
- Safe for patient and user, thanks to the adjustable sliding depth stopper
- Graduated puncture cannula to assist in precise depth placement of the needle
- Transparent luer cone for visual inspection of aspirate
- Colour coded hub for easy identification of needle gauge size



Reusable syringe gun (code: SG01)

### CYTOCAN Qty/Box 20pcs

CODE	SIZE
● VVB1810	18Gx10cm
● VVB1815	18Gx15cm
● VVB1820	18Gx20cm
● VVB2010	20Gx10cm
● VVB2015	20Gx15cm
● VVB2020	20Gx20cm
● VVB2105	21Gx5cm
● VVB2110	21Gx10cm
● VVB2115	21Gx15cm
● VVB2120	21Gx20cm
● VVB2210	22Gx10cm
● VVB2215	22Gx15cm
● VVB2220	22Gx20cm
● VVB2310	23Gx10cm
● VVB2315	23Gx15cm
● VVB2505	25Gx5cm

On request Vigeo is available to provide needles with different diameter and length





- Individually packaged in sterile wrappers
- Intended for use in amniocentesis or chorionic villus sampling
- Safe for patient and user, thanks to the adjustable sliding depth stopper
- Graduated puncture cannula to assist in precise depth placement of the needle
- Transparent luer cone for visual inspection of aspirate
- Colour coded hub for easy identification of needle gauge size
- Sharp, non-coring needle tip to minimize patient trauma
- Self-locating stylets to prevent blockage

## AMNIO



## CHORIAL KIT



## PRENATAL KIT



### AMNIO Qty/Box 20pcs

CODE	SIZE
● VVBP1915	19Gx15cm
● VVBP1920	19Gx20cm
● VVBP2015	20Gx15cm
● VVBP2020	20Gx20cm
● VVBP2115	21Gx15cm
● VVBP2120	21Gx20cm
● VVBP2315	23Gx15cm

### CHORIAL KIT Qty/Box 20pcs

● VVBU1815	18Gx15cm
------------	----------

### PRENATAL KIT Qty/Box 20pcs

● VVBL2016	20Gx16cm
● VVBL2020	20Gx20cm
● VVBL2116	21Gx16cm
● VVBL2120	21Gx20cm

On request Vigeo is available to provide needles with different diameter and length



- Designed for quick and easy cortical penetration
- Ergonomic T-handle for increased control
- Stylet with pyramidal tip
- Cannula with fish mouth tip or with crown tip
- Centimeter depth markings assist in needle placement
- Spoon collection cannula or forceps collection cannula
- Improved patient comfort



VVJTxyy  
Spoon collection cannula



VVJTxyy10  
Forceps collection cannula

### BLOCK Qty/Box 10pcs

**VVJTxyy** Fish mouth tip and spoon collection cannula

**VVJTxyy10** Crown tip and forceps collection cannula

CODE	SIZE
● VVJT0810 VVJT081010	8Gx10cm
● VVJT0815 VVJT081510	8Gx15cm
● VVJT0910 VVJT091010	9Gx10cm
● VVJT0915 VVJT091510	9Gx15cm
● VVJT1110 VVJT111010	11Gx10cm
● VVJT1115 VVJT111510	11Gx15cm
● VVJT1307 VVJT130710	13Gx7cm
● VVJT1310 VVJT131010	13Gx10cm

On request Vigeo is available to provide needles with different diameter and length



- Designed for quick and easy cortical penetration
- Ergonomic T-handle for increased control
- Stylet with pyramidal tip
- Cannula with fish mouth tip or with crown tip
- Centimeter depth markings assist in needle placement
- Extraction cannula with centimeter markings and blunt tip
- Improved patient comfort



Side distal holes (VVJExxy)



Fish mouth tip

### BCUT Qty/Box 10pcs

**VVJxxyy** Fish mouth tip collection cannula

**VVJxxyy10** Side distal hole for bone marrow aspiration

**VVJxxyy20** Crown tip cannula

**VVJExxyy** Side distal holes for efficient aspiration during harvesting procedures

CODE				SIZE
●	VVJ0810	VVJ081010	VVJ081020	8Gx10cm
●	VVJ0815	VVJ081510	VVJ081520	8Gx15cm
●	VVJ0910	VVJ091010	VVJ091020	9Gx10cm
●	VVJ0915	VVJ091510	VVJ091520	9Gx15cm
●	VVJ1110	VVJ111010	VVJ111020	VVJE1110
●	VVJ1115	VVJ111510	VVJ111520	VVJE1115
●	VVJ1307	VVJ130710	VVJ130720	13Gx7cm
●	VVJ1310	VVJ131010	VVJ131020	VVJE1310
				13Gx10cm

Bone and bone marrow biopsy device



- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



Side distal holes VVIE1xxyy/xx(x)

**MYELO SLIM** Qty/Box 10pcs

**VVIE1xxyy/xx(x)** for bone marrow harvesting needle

CODE		SIZE
VVI11410/50	VVIE11410/50	14Gx30mm
VVI11430/70	VVIE11430/70	14Gx50mm
VVI11510/50	VVIE11510/50	15Gx30mm
VVI11530/70	VVIE11530/70	15Gx50mm
VVI11610/50	VVIE11610/50	16Gx30mm
VVI11630/70	VVIE11630/70	16Gx50mm
VVI11810/50	VVIE11810/50	18Gx30mm

On request Vigeo is available to provide needles with different diameter and length

Bone and bone marrow biopsy device



- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



Side distal holes VVIExxxy/xx(x)

MYELO Qty/Box 10pcs		
VVIExxxy/xx(x) for bone marrow harvesting needle		
CODE		SIZE
VVI1410/50	VVIE1410/50	14Gx30mm
VVI1430/70	VVIE1430/70	14Gx50mm
VVI1510/50	VVIE1510/50	15Gx30mm
VVI1530/70	VVIE1530/70	15Gx50mm
VVI1610/50	VVIE1610/50	16Gx30mm
VVI1630/70	VVIE1630/70	16Gx50mm
VVI1810/50	VVIE1810/50	18Gx30mm

On request Vigeo is available to provide needles with different diameter and length


Bone and bone marrow biopsy device



- Ergonomic T-handle design to minimize hand and arm fatigue
- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



Side distal holes VVILExxyy/xx(x)



Colour code



ERGO Qty/Box 10pcs		
VVILExxyy/xx(x) for bone marrow harvesting needle		
CODE		SIZE
VVIL1410/50	VVILE1410/50	14Gx30mm
VVIL1430/70	VVILE1430/70	14Gx50mm
VVIL1510/50	VVILE1510/50	15Gx30mm
VVIL1530/70	VVILE1530/70	15Gx50mm
VVIL1610/50	VVILE1610/50	16Gx30mm
VVIL1810/50	VVILE1810/50	18Gx30mm

On request Vigeo is available to provide needles with different diameter and length



# D PUNCH/CURETTE SKIN CARE

CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )

Disposable skin biopsy device



- Sterile, single-use, stainless steel instrument for accurate and gentle skin biopsies
- Stainless steel sharp cutting edge
- Ergonomic ribbed handle and anti-slip surface for sure grip control
- Size embossed on each punch for easy identification
- Wide range of sizes



Ergonomic ribbed handle



Curette version (VPTC)

## D PUNCH/CURETTE Qty/Box 20 pcs

CODE	SIZE
VPT2 VPTC2	2,0mm
VPT3 VPTC3	3,0mm
VPT4 VPTC4	4,0mm
VPT5 VPTC5	5,0mm
VPT6 —	6,0mm
VPT7 VPTC7	7,0mm
VPT8 —	8,0mm

On request Vigeo is available to provide needles with different diameter and length

Multilayer decontamination mat















- Used in all environments that require the total absence of contaminating agents and the preservation of high hygiene levels
- Barrier against the growth and spreading of microorganisms
- Antibacterial and antifungal glue captures and retains dust and bacteria
- 1,4mm thickness mat includes 30 adhesive PE sheets
- Non-adhesive strip for easy removal of used layers
- Dispenser folder for comfortable handling and vertical storage




#### DECOMAT Qty/Box 5pcs

CODE	SIZE
● DK-1245-BFVG	115x45 cm
● DK-1260-BFVG	115x60 cm
● DK-1290-BFVG	115x90 cm
● DK-1245-VFVG	115x45 cm
● DK-1260-VFVG	115x60 cm
● DK-1290-VFVG	115x90 cm

On request Vigeo is available to provide different measures

14 G	=		2,10 mm
16 G	=		1,65 mm
17 G	=		1,47 mm
18 G	=		1,27 mm
19 G	=		1,10 mm
20 G	=		0,91 mm
21 G	=		0,82 mm
22 G	=		0,72 mm
23 G	=		0,64 mm
24 G	=		0,55 mm
25 G	=		0,51 mm
27 G	=		0,39 mm





Urology

Interventional

Tissue Marker

Biopsy Solutions

Radiology

Breast care

R&D

**VIGEO** srl  
 Via Dell' Alpino 2  
 46031 Bagnolo San Vito MN - ITALY  
 Tel. +39 0376 449097 - Fax +39 0376 282245  
 export@vigeosrl.it

ISO 13485:2016  
[www.vigeosrl.it](http://www.vigeosrl.it)  
 September 2023

